

MICROPURE Monofocal Hydrophobic IOL



SEE FURTHER FOR YOUR CATARACT PATIENTS



MICROPURE & MICROPURE 123

See further for your cataract patients

The MICROPURE lens (MICAGF (123I) is a monofocal intraocular implant based on proven 4 closed loops platform stability with a reliable hydrophobic acrylic raw material.

The Micropure is also available preloaded in a cartridge, which is clipped to the Single-Use injector 1.2.3 Premium.



MICROPURE 123

Where Micro refers to "support Corresponds with the four closed loops haptics"

hydrophobic material (GFY)

Refers to the preloaded system




BASED ON MICRO HAPTIC DESIGN

The closed quadripode Micro haptic design has a widely adapted geometry, which dates back to 2003.

Besides, the large contact angle of the 4 haptics supports the lens centration at different simulated capsular bag diameters in vitro.1

Ø capsule: 10,0mm	Micropure	1-piece C-loop yellow lens	1-piece C-loop lens
Contact Angle	135°	61°	57°



Optimized stability thanks to 4 points of contact

Under the capsular retraction effect, the haptics absorb the forces applied on the periphery to maintain the implant optic correctly centred²

Compressed implant \emptyset 9.50mm




1 ISO 11979-3 VALIDATION REPORT: MECHANICAL PROPERTIES – Report PhysIOL October 2015 | 2 Mechanical properties according to ISO 11979-3 — PhysIOL test report dated 23 June 2014 | 3 Glistenings in Alcon Acrysof Intraocular Lenses, GEORGE H.H. BEIKO, B.M., B.Ch., FRCSCST. CATHARINES, CANADAASSIST PROF, MCMASTER UNIVLECTURER, UNIV OF TORONTO (https://slideplayer.com/slide/10336530)| 4 Miyata A, Jpn J Ophthalmol 2001, 45(6):564-569.| 5 https://www.eurotimes.org/capsular-bag-stability-findl| 6 CER F2 (MIC-GFY) | RD-REP-210-1-2021 | V1.0 | 27.04.2021 | 7 Biomaterial Optical Purity. The David J Apple International Laboratory for Ocular Pathology, 3 MAY 2017.


"Why implant a lens when the potential for glistening exists, if there is the equivalent-quality IOL available that carries out all the function of the lens but is not encumbered with glistening." - David Apple cited by Georges Beiko³

BASED ON HYDROPHOBIC (GFY)


Proven 10 years reliable clinical outcomes from unique G-free hydrophobic

The GFY raw material is a Grade 0 based on the Miyata glistening scale (in vitro) 4, 5, 6





3



Contrast Sensitivity [llogCS] 2.5 2.0 1.5 1.0 3M, N=18 _ _ population norm 0.5 0.0

12

Spatial Frequency [cpd]

18

3

BASED ON CLINICAL EVALUATION^{6,7}

- 1 Pre- and postoperative MICROPURE 123 IOLs)
- 2 At 3 months, monocular photopic DCIVA was 0.50 logMAR (20/63) or better in 42.9% and 33.3% of the eyes respectively.
 - 3 Mean binocular photopic contrast sensitivity at 3 months.



PhysIOL MICROPURE 123

Monofocal Hydrophobic Preloaded



Technical Specifications

Commercial name	MICROPURE 123				
Material	PhysIOL G-free® (GFY) (hydrophobic acrylic glistening-free)¹				
Overall diameter	0D to 24.5D: 11.00 mm 25D to 30D: 10.75 mm				
Optic diameter	0D to 24.5D: 6.00 mm 25D to 30D: 5.75 mm				
Optic	Aspheric aberration-correcting (-0.11µ SA)				
Filtration	UV & blue light				
Refractive index	1.52				
Abbe number	42				
Angulation	2°				
Injection system	PhysIOL 1.2.3				
Incision size	≥ 2.2 mm				
Spherical power	0D to 9D (1D steps) & 10D to 30D (0.5D steps) Cartridge with PRS® technology²				
Square edge	360°				
Nominal manufacturer A constant	119.40				
Suggested A constant ³		Interferometry	Ultrasound		
	Hoffer Q: pACD	5.85	5.59		
	Holladay 1: Sf	2.06	1.80		
	Barrett: LF	2.09	-		
	SRK/T: A	119.40	119.05		
	Haigis⁴: a0; a1; a2	1.70; 0.4; 0.1	1.214; 0.4; 0.1		
	MICROPURE (non-preloaded)				
Overall diameter	-10D to 24.5D: 11.00 mm & 25D to 35D: 10.75 mm				
Optic diameter	-10D to 24.5D: 6.00 mm & 25D to 35D: 5.75 mm				
Injection system	Medicel Accuject 1.8 up to 24.5D & Accuject 2.0/2.1/2.2 up to 35D				
Incision size	≥ 1.8 mm				
Spherical power	-10D to 9D (1D steps) & 10D to 30D (0.5D steps) & 31D to 35D (1D steps)				

³ Estimates only: surgeons are recommended to use their own values based upon their personal experience. Refer to our website for updates. ⁴ Not optimized.





¹ The PhysIOL G-free® (GFY) is patented since 2010. Chassain C, J Fr Ophthalmol 2018, 41(6):513-520.

 $^{^2\,\}mbox{The PRS}^{\tiny \textcircled{\tiny 0}}$ technology is patent pending.